

INSTRUCTION FOR USE

DISPOSAFE Safety IV Cannula with wings injection port

1. MATERIALS USED:

2. COATING MATERIAL: Polydimethylsiloxane

3. DESCRIPTION: Each IV Cannula with safety features is a device that can be used for infusion/ administration through peripheral vein. The device is enabled with passive safety feature, to eliminate the needle-stick injuries to the healthcare professionals and waste handlers. Inherently designed transparent flashback chamber enable ease of conformation to venous access. The device is enable with female 6% taper, facilitates integrity of connection with other combination devices.

4. INDICATIONS:

These instructions contain important information for use of the product. Read the entire contents of these instructions for use, including warnings, cautions and falling to follow instructions could result in death or serious injury to the patient and/or clinician.

PP, PE, POM, HDPE, ABS, Silicon Rubber, Stainless Steel, FEP/PURIPTFE.

4.1. The device is intended to be used for sample blood or administer fluids intravenously.

4.2. The device may be used on any patient population with consideration given to adequacy of vascular anatomy and appropriateness of procedure.

4.3. Intermittent intravenous drug administration.

4.4. Prophylactic creation of a secure venous access in patients who may require urgent intravenous drug administration, in particular prior to diagnostic or therapeutic procedures

5. CONTRAINDICATIONS:

5.1. The patient is known or is suspected to be allergic to materials contained in the device.

5.2. Past irradiation of prospective insertion site.

5.3. Previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.

6. INSTRUCTION FOR USE:

6.1. Carefully select and aseptically prepare the site.

6.2. Select suitable size of I.V. Cannula & inspect visually to ascertain that package is intact.

6.3. Remove the protective sheath in straight outward motion and inspect the device. Ensure catheter hub

6.4. Approach the insertion site slowly at low angle.

6.5. Confirm successful venipuncture by visualizing blood in the flashback chamber.

6.6. Advance the catheter further into the vein, such that catheter is nearly parallel to skin surface, while slightly withdrawing the steel needle.

6.7. Once the insertion is complete, press the vein above from the catheter tip, and pull the needle straight outward. Safety clip will automatically attach to needle tip as needle tip exits catheter hub. Confirm the activation of safety mechanism by hearing an audible click.

6.8. Hold the catheter body and detach the safety cage.

6.9. Connect the infusion line and cover the puncture site with sterile dressing.

6.10. Fix an appropriate sterile lock to avoid the contamination and spillage of blood, while not in use or

6.11. Bolus injection is possible via the integrated injection valve."

7. CAUTION:

7.1. adhered to avoid the exposure to blood borne pathogens. It is recommended that the healthcare professionals follow the recommendations set forth by the Centers for Disease Control and

. The device is designed to reduce the risk of accidental needle-sticks; however, care must be taken in order to avoid needle-sticks, in addition to the instruction for use. Universal precautions must be

Prevention/Occupational Safety and Health Administration (CDC/OSHA).

7.2. The catheter is soft and fragile; do not use any sharp equipment near to the insertion site or during the insertion or removal.

8. WARNINGS:

8.1. In order to prevent misuse of the device, the use of this product is restricted to a Doctor or qualified medical practitioner.

8.2. The product should be used according to the instructions for use.

8.3. If there is any change in expected performance of the device or in case of any malfunction the device.

8.4. Any device that is connected to this product must comply with ISO 80369-7, in order to achieve the

intended performance of this product & to avoid leakage in the connection.

8.5. The product should not be reprocessed.

8.6. Visually inspect and carefully check the product and packaging before use. Improper transport and handling may

cause structural and/or functional damage to device or packaging

8.7. The product is guaranteed sterile & non-pyrogenic, if the package has not been opened or damaged.

8.8. Do not clean or resterilize. For single use only. Discard after use.

8.9. Do not expose to heat or direct sunlight

8.10. The product should be used immediately after opening the packaging.

8.11. Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections.

















8.12. Disposal Discard: Dispose of/ Discard the used Device in accordance with your Country's Healthcare and Safety Regulations.


8.13. When using alcohol or alcohol containing antiseptics, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing.

8.14. Alcohol should not be used to lock, soak, or decontaminate the catheters because alcohol is known to degrade the catheter properties over time with repeated and prolonged exposure.


8.15. (Pediatric) Insertion

techniques and placement locations are often modified according to the size and

	Do Not Reuse		Lot Number
	Do Not Resterilize		Date of Manufacturing
	Do not use if package is damaged and consult instruction for use.		Use By / Expiry Date
	Reference No.		Consult Instructions For Use
	Manufacturer		Keep dry
	Sterilized by Ethylene oxide Single Sterile Barrier System		Keep away from sunlight
	Non-Pyrogenic		Medical Device
	Temperature Limit		Caution




1. Remove The Cannula Cap




2. Hold The Bevel Phase Of The Cannula Upwards.




3. Insert Cannula Into Vein At A Low Angle.




4. Withdraw The Needle Gently.



5. Safety Clip Covered Needle Tip.



6. Attach The Threaded Stopper



7. Dispose of the needle in sharps container.

EC REP

CEPARTNER4U B.V., Edoomlaan 13,
3951 DB Maar, The Netherlands.
E-mail : Office@cepartner4u.com
AW/IFU_032, REV:09

Disposafe Health and Life Care Limited

Plot No.1 & 2, Phase-II, Sector-59, Ballabgarh,

Faridabad -121004, Haryana, INDIA

E-mail:info@disposafe.com

Web:www.disposafe.com

Customer Care No. +91-129-4170201

Mfg. Lic. No. MFG/MD/2019/000106

CE
0123